SORIMEX

Kable czujników temperatury
Cable for temperature sensors

SORIMEX sp. z o.o. sp. k.
ul. Równinna 25, 87-100 Toruń, Poland
tel.: +48 56 657 77 20, fax: +48 56 657 77 21
e-mail: biuro@sorimex.pl ◆ www.sorimex.eu

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EN

Instruction of use and storage cable for temperature sensor

1. General description
The temperature sensor cable is intended for the electrical connection of the temperature sensor to the device. The temperature sensor cable is classified as class I and is supplied as non-sterile.

2. Characteristics
The temperature sensor cable consists of a suitable multi-core cable and its connector is adapted to the temperature sensor and the device to which it is to cooperate. The user should make sure that the cable is suitable for the specific type of temperature sensor and device .
Cables are intended for use with undamaged patient's skin.

3. Principles of use
When using cables for temperature sensor, the following principles should be observed:
- make sure that cable is suitable for the device to work with,
- cables should not be excessively bent (bend radius should not be lower than 60 mm),
- do not pull or jerk the cable,
- place the cable so that there would be no risk of mechanical damage, e.g. by stepping on, running over, placing a heavy/sharp object on it, cutting etc.,
- hold the connector cover (socket, pin) when connecting and disconnecting the cable,
- protect the cable connectors against suffusion.

4. Installation
Connect the temperature sensor cable from one side to the temperature sensor and the other side to the device. Turn on the unit and check that the monitor operates properly according to the monitor's instruction manual.Pay attention to the proper selection of the temperature sensor cable for the appropriate device and temperature sensor.

5. Recommended cleaning, disinfection and sterilisation
These activities should be realized every time before use of the device and cable before the examination of the next patient. Before starting the cleaning, disinfection, or sterilisation, the cable should be disconnected from the device with which it cooperates.

Cleaning
Wash the cable with warm soap water, delicately clean using soft cloth wetted with soap water (mild detergent) or 70% isopropyl alcohol. The cable should dry before use.
Do not use corrosive substances or other chemical preparations.
When cleaning, do not soak or immerse the cables, do not wash under running water, and do not use excessive strength which may rupture the internal conductors, leading to cable damage. Continuous bending of the conductors during use and cleaning may rupture the internal conductors. Damage of this type is not covered by the guarantee.
Avoid contact with strong solvents: aromatic, chloric, ketone, etheral, and ester. Do not allow for contact of cleaners with tplugs and pins on either side of the cable.

Disinfection
Cables should be disinfected by wiping with 70% isopropyl alcohol and leaving them to dry. Do not allow for contact of cleaners with tplugs and pins on either side of the cable.

Sterilization
The recommended method consists of sterilization using ethylene oxide. After sterilization and before use, cables should be accurately degassed in device intended for that purpose.
Sterilization with other methods is not recommended for it may permanently damage the cable.

6. Packaging and storing
Temperature sensor cables are packed individually. In order to maximize the product's working life, they should be stored in the original packaging, in accordance with provided storing conditions.

Storage conditions:
- temperature: -25°C to + 55°C,
- protect against sunlight, moisture and dirt.

Cables in use should be stored in a dry place.
When not used, temperature sensor cables should be loosely coiled, without sharp bends.

7. Warnings ⚠
Failure to observe the above rules may result in a break in the electrical circuit in the transmission cable or cause a short circuit. These symptoms are easy to identify by operating the camera. The damaged cable must be replaced and replaced with a working one.

Utilisation of used cable should be done in accordance with binding provisions.

8. Guarantee
The guarantee period covers 1 year from the date of sale.
The guarantee does not include defects resulting from improper use and storage.

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